

## Copper Setting-up Sample, 50 g Chips

Art. ID MBH-165XALB1SUSN-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	0,03	%		
	Lead (Pb)	[7439-92-1]	0,2	%		
	Zinc (Zn)	[7440-66-6]	0,06	%		
	Iron (Fe)	[7439-89-6]	2,8	%		
	Nickel (Ni)	[7440-02-0]	5,3	%		
	Aluminium (Al)	[7429-90-5]	9	%		
	Silicon (Si)	[7440-21-3]	0,1	%		
	Arsenic (As)	[7440-38-2]	0,005	%		
	Manganese (Mn)	[7439-96-5]	0,08	%		
	Bismuth (Bi)	[7440-69-9]	0,015	%		
	Phosphorus (P)	[7723-14-0]	0,015	%		
	Chromium (Cr)	[7440-47-3]	0,01	%		
	Magnesium (Mg)	[7439-95-4]	0,04	%		
	Copper (Cu)	[7440-50-8]	82	%		