

Inconel 600 - Ni/Cr/Fe Alloy, 50 g Chips

Art. ID MBH-28X6003E-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,025	%		
	Cobalt (Co)	[7440-48-4]	0,62	%		
	Chromium (Cr)	[7440-47-3]	15,56	%		
	Copper (Cu)	[7440-50-8]	0,42	%		
	Iron (Fe)	[7439-89-6]	7,1	%		
	Magnesium (Mg)	[7439-95-4]	0,01	%		
	Manganese (Mn)	[7439-96-5]	0,47	%		
	Silicon (Si)	[7440-21-3]	0,74	%		
	Titanium (Ti)	[7440-32-6]	0,22	%		