

Leaded Brass Alloy, 50 g Chips

Art. ID MBH-31X7835.1U-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	0,446 ± 0,005	%		
	Lead (Pb)	[7439-92-1]	2,78 ± 0,02	%		
	Zinc (Zn)	[7440-66-6]	33,42 ± 0,05	%		
	Iron (Fe)	[7439-89-6]	0,173 ± 0,005	%		
	Nickel (Ni)	[7440-02-0]	0,12 ± 0,003	%		
	Aluminium (Al)	[7429-90-5]	0,0187 ± 0,0009	%		
	Manganese (Mn)	[7439-96-5]	0,0105 ± 0,0005	%		
	Silicon (Si)	[7440-21-3]	0,0218 ± 0,0014	%		
	Phosphorus (P)	[7723-14-0]	0,0493 ± 0,0013	%		
	Cobalt (Co)	[7440-48-4]	0,0091 ± 0,0003	%		
	Bismuth (Bi)	[7440-69-9]	0,0089 ± 0,0008	%		
	Antimony (Sb)	[7440-36-0]	0,0156 ± 0,0009	%		
	Arsenic (As)	[7440-38-2]	0,0138 ± 0,0009	%		
	Chromium (Cr)	[7440-47-3]	0,0028 ± 0,0002	%		
	Cadmium (Cd)	[7440-43-9]	0,0037 ± 0,0003	%		
	Silver (Ag)	[7440-22-4]	0,0083 ± 0,0004	%		
	Tellurium (Te)	[13494-80-9]	0,0015 ± 0,0001	%		
	Copper (Cu)	[7440-50-8]	62,9 ± 0,08	%		