

Leaded Brass Alloy, 50 g Chips

Art. ID MBH-31X7835.8B-C
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Lead (Pb) | [7439-92-1] | 3,22 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,451 | % | | |
| | Zinc (Zn) | [7440-66-6] | 21,55 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,0446 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,157 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,0102 | % | | |
| | Aluminium (Al) | [7429-90-5] | 0,219 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,151 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,101 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,11 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,154 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,313 | % | | |
| | Silver (Ag) | [7440-22-4] | 0,549 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,0944 | % | | |
| | Copper (Cu) | [7440-50-8] | 72,7 | % | | |