

## **Brass - Cu/Zn Binary, 50 g Chips**

Art. ID MBH-31XB3N-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0152 ± 0,0006	%		
	Arsenic (As)	[7440-38-2]	0,0104 ± 0,0007	%		
	Boron (B)	[7440-42-8]	0,0008 ± 0,0001	%		
	Bismuth (Bi)	[7440-69-9]	0,0152 ± 0,0005	%		
	Cadmium (Cd)	[7440-43-9]	0,0027 ± 0,0001	%		
	Cobalt (Co)	[7440-48-4]	0,0066 ± 0,0004	%		
	Copper (Cu)	[7440-50-8]	64,56 ± 0,09	%		
	Iron (Fe)	[7439-89-6]	0,0408 ± 0,0012	%		
	Manganese (Mn)	[7439-96-5]	0,029 ± 0,0008	%		
	Nickel (Ni)	[7440-02-0]	0,0171 ± 0,0004	%		
	Phosphorus (P)	[7723-14-0]	0,0366 ± 0,0013	%		
	Lead (Pb)	[7439-92-1]	0,085 ± 0,003	%		
	Antimony (Sb)	[7440-36-0]	0,0148 ± 0,0007	%		
	Tin (Sn)	[7440-31-5]	0,0394 ± 0,0012	%		
	Tellurium (Te)	[13494-80-9]	0,0056 ± 0,0005	%		
	Zinc (Zn)	[7440-66-6]	35,05 ± 0,1	%		