

High-Tensile Brass Alloy, 50 g Chips

Art. ID MBH-31XHT31B-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	6,82	%		
	Arsenic (As)	[7440-38-2]	0,0005	%		
	Boron (B)	[7440-42-8]	0,0014	%		
	Carbon (C)	[7440-44-0]	0,0057	%		
	Copper (Cu)	[7440-50-8]	67	%		
	Iron (Fe)	[7439-89-6]	3,01	%		
	Manganese (Mn)	[7439-96-5]	5,69	%		
	Nickel (Ni)	[7440-02-0]	0,226	%		
	Phosphorus (P)	[7723-14-0]	0,003	%		
	Lead (Pb)	[7439-92-1]	0,0139	%		
	Sulfur (S)	[7704-34-9]	0,0007	%		
	Silicon (Si)	[7440-21-3]	0,0443	%		
	Tin (Sn)	[7440-31-5]	0,079	%		
	Zinc (Zn)	[7440-66-6]	17,06	%		