

Phosphor Bronze Alloy, 50 g Chips

Art. ID MBH-32X52100A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0011	%		
	Aluminium (Al)	[7429-90-5]	0,0009	%		
	Arsenic (As)	[7440-38-2]	0,0009	%		
	Bismuth (Bi)	[7440-69-9]	0,0019	%		
	Copper (Cu)	[7440-50-8]	92,1	%		
	Iron (Fe)	[7439-89-6]	0,0008	%		
	Nickel (Ni)	[7440-02-0]	0,0111	%		
	Phosphorus (P)	[7723-14-0]	0,146	%		
	Lead (Pb)	[7439-92-1]	0,0031	%		
	S (Sulfur)		0,0008	%		
	Tin (Sn)	[7440-31-5]	7,73	%		
	Zinc (Zn)	[7440-66-6]	0,0026	%		