

## Gun Metal Alloy, 50 g Chips

Art. ID MBH-33XGM20B-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Tin (Sn)          | [7440-31-5] | 4,45                | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,29                | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 1,78                | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,43                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,21                | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 0,13                | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,29                | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,04                | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,04                | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 2,35                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,06                | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 0,02                | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,02                | %    |        |        |
|                  | Silver (Ag)       | [7440-22-4] | 0,21                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | [bal]               | %    |        |        |