

## Gun Metal Alloy, 50 g Chips

Art. ID MBH-33XGM4AD-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	3,02	%		
	Lead (Pb)	[7439-92-1]	5,27	%		
	Zinc (Zn)	[7440-66-6]	5,9	%		
	Iron (Fe)	[7439-89-6]	0,0932	%		
	Nickel (Ni)	[7440-02-0]	1,482	%		
	Aluminium (Al)	[7429-90-5]	0,0015	%		
	Silicon (Si)	[7440-21-3]	0,001	%		
	Arsenic (As)	[7440-38-2]	0,0228	%		
	Manganese (Mn)	[7439-96-5]	0,0006	%		
	Bismuth (Bi)	[7440-69-9]	0,0442	%		
	Antimony (Sb)	[7440-36-0]	0,0568	%		
	Phosphorus (P)	[7723-14-0]	0,0034	%		
	Sulfur (S)	[7704-34-9]	0,034	%		
	Cobalt (Co)	[7440-48-4]	0,0077	%		
	Silver (Ag)	[7440-22-4]	0,0206	%		
	Copper (Cu)	[7440-50-8]	84,02	%		