

## Gun Metal Alloy, 50 g Chips

Art. ID MBH-33XGM9A-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Silver (Ag)       | [7440-22-4] | 0,025               | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 0,005               | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,025               | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,075               | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,07                | %    |        |        |
|                  | Cobalt (Co)       | [7440-48-4] | 0,08                | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | Rem                 | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,09                | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,73                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,055               | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 7                   | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,18                | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,002               | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 3                   | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | 13,7                | %    |        |        |