

### **Cu/Ni/Sn (Spinodal Alloy), 50 g Chips**

Art. ID MBH-36XSP1A-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,005	%		
	Aluminium (Al)	[7429-90-5]	0,002	%		
	Boron (B)	[7440-42-8]	0,0007	%		
	Bismuth (Bi)	[7440-69-9]	0,0039	%		
	Cobalt (Co)	[7440-48-4]	0,057	%		
	Copper (Cu)	[7440-50-8]	84,9	%		
	Iron (Fe)	[7439-89-6]	0,45	%		
	Manganese (Mn)	[7439-96-5]	0,084	%		
	Nickel (Ni)	[7440-02-0]	8,33	%		
	Lead (Pb)	[7439-92-1]	0,0115	%		
	Sulfur (S)	[7704-34-9]	0,005	%		
	Antimony (Sb)	[7440-36-0]	0,0177	%		
	Silicon (Si)	[7440-21-3]	0,004	%		
	Tin (Sn)	[7440-31-5]	5,75	%		
	Zinc (Zn)	[7440-66-6]	0,344	%		