

Cu/Ni/Sn (Spinodal Alloy), 50 g Chips

Art. ID MBH-36XSP2A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0181	%		
	Aluminium (Al)	[7429-90-5]	0,0003	%		
	Boron (B)	[7440-42-8]	0,0005	%		
	Cobalt (Co)	[7440-48-4]	0,119	%		
	Copper (Cu)	[7440-50-8]	74,91	%		
	Magnesium (Mg)	[7439-95-4]	0,0002	%		
	Manganese (Mn)	[7439-96-5]	0,0019	%		
	Nickel (Ni)	[7440-02-0]	15,72	%		
	Lead (Pb)	[7439-92-1]	0,026	%		
	Sulfur (S)	[7704-34-9]	0,003	%		
	Antimony (Sb)	[7440-36-0]	0,006	%		
	Tin (Sn)	[7440-31-5]	8,92	%		
	Zinc (Zn)	[7440-66-6]	0,029	%		