

## **Al/Mg Alloy, 50 g Chips**

Art. ID MBH-511XG3000B1C-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	(<0,0001)	%		
	Chromium (Cr)	[7440-47-3]	0,0972	%		
	Copper (Cu)	[7440-50-8]	0,287	%		
	Iron (Fe)	[7439-89-6]	0,669	%		
	Gallium (Ga)	[7440-55-3]	0,051	%		
	Magnesium (Mg)	[7439-95-4]	0,253	%		
	Manganese (Mn)	[7439-96-5]	0,752	%		
	Nickel (Ni)	[7440-02-0]	0,142	%		
	Lead (Pb)	[7439-92-1]	0,145	%		
	Silicon (Si)	[7440-21-3]	0,788	%		
	Tin (Sn)	[7440-31-5]	0,163	%		
	Titanium (Ti)	[7440-32-6]	0,0327	%		
	Vanadium (V)	[7440-62-2]	0,054	%		
	Zinc (Zn)	[7440-66-6]	0,0347	%		