

Al/Mg Alloy with Bi, 50 g Chips

Art. ID MBH-511XG6065A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	0,0011	%		
	Bismuth (Bi)	[7440-69-9]	1,27	%		
	Chromium (Cr)	[7440-47-3]	0,0759	%		
	Copper (Cu)	[7440-50-8]	0,214	%		
	Iron (Fe)	[7439-89-6]	0,74	%		
	Magnesium (Mg)	[7439-95-4]	1,18	%		
	Manganese (Mn)	[7439-96-5]	0,131	%		
	Nickel (Ni)	[7440-02-0]	0,0144	%		
	Lead (Pb)	[7439-92-1]	0,031	%		
	Silicon (Si)	[7440-21-3]	0,495	%		
	Tin (Sn)	[7440-31-5]	0,0251	%		
	Titanium (Ti)	[7440-32-6]	0,112	%		
	Zinc (Zn)	[7440-66-6]	0,164	%		
	Zirconium (Zr)	[7440-67-7]	0,144	%		