

Mg/Al Alloy, 50 g Chips

Art. ID MBH-64XMgQ1A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1,083	%		
	Beryllium (Be)	[7440-41-7]	0,00036	%		
	Copper (Cu)	[7440-50-8]	0,084	%		
	Iron (Fe)	[7439-89-6]	0,0034	%		
	Manganese (Mn)	[7439-96-5]	0,377	%		
	Nickel (Ni)	[7440-02-0]	0,0041	%		
	Lead (Pb)	[7439-92-1]	0,02	%		
	Silicon (Si)	[7440-21-3]	0,062	%		
	Tin (Sn)	[7440-31-5]	0,0195	%		
	Zinc (Zn)	[7440-66-6]	0,235	%		