

Tin-base (Lead-Free) Solder, 50 g Chips

Art. ID MBH-74XBZ1A-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,004	%		
	Aluminium (Al)	[7429-90-5]	0,0021	%		
	Arsenic (As)	[7440-38-2]	0,0119	%		
	Bismuth (Bi)	[7440-69-9]	3,03	%		
	Cadmium (Cd)	[7440-43-9]	0,0012	%		
	Copper (Cu)	[7440-50-8]	0,026	%		
	Iron (Fe)	[7439-89-6]	0,011	%		
	Nickel (Ni)	[7440-02-0]	0,0097	%		
	Lead (Pb)	[7439-92-1]	0,0238	%		
	Antimony (Sb)	[7440-36-0]	0,031	%		
	Zinc (Zn)	[7440-66-6]	8,27	%		