

***** sold out, please see MBH-74XCA4-21 *** Tin-base (Lead-Free) Solder, 50 g Chips**

Art. ID MBH-74XCA4D-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	3,07 ± 0,05	%		
	Arsenic (As)	[7440-38-2]	0,003 ± 0,0003	%		
	Gold (Au)	[7440-57-5]	0,002 ± 0,0001	%		
	Bismuth (Bi)	[7440-69-9]	0,0515 ± 0,0011	%		
	Cadmium (Cd)	[7440-43-9]	0,0025 ± 0,0001	%		
	Cobalt (Co)	[7440-48-4]	0,0005 ± 0,0001	%		
	Chromium (Cr)	[7440-47-3]	0,001 ± 0,0002	%		
	Copper (Cu)	[7440-50-8]	0,552 ± 0,007	%		
	Iron (Fe)	[7439-89-6]	0,004 ± 0,0003	%		
	Mercury (Hg)	[7439-97-6]	0,0075 ± 0,0006	%		
	Indium (In)	[7440-74-6]	0,009 ± 0,001	%		
	Nickel (Ni)	[7440-02-0]	0,0589 ± 0,0012	%		
	Lead (Pb)	[7439-92-1]	0,081 ± 0,003	%		
	Antimony (Sb)	[7440-36-0]	0,0612 ± 0,0017	%		
	Zinc (Zn)	[7440-66-6]	0,0051 ± 0,0005	%		