

## **Tin-base (Lead-Free) Solder, 50 g Chips**

Art. ID MBH-74XHAG-C  
Unit 50 g  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	2,8	%		
	Arsenic (As)	[7440-38-2]	0,0032	%		
	Bismuth (Bi)	[7440-69-9]	0,0639	%		
	Cadmium (Cd)	[7440-43-9]	0,0018	%		
	Copper (Cu)	[7440-50-8]	0,629	%		
	Iron (Fe)	[7439-89-6]	0,0029	%		
	Indium (In)	[7440-74-6]	0,009	%		
	Nickel (Ni)	[7440-02-0]	0,0133	%		
	Lead (Pb)	[7439-92-1]	0,077	%		
	Antimony (Sb)	[7440-36-0]	2,1	%		
	Zinc (Zn)	[7440-66-6]	2,73	%		