

Tin-base (Lead-Free) Solder, 50 g Chips

Art. ID MBH-74XHBG-C
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,086	%		
	Aluminium (Al)	[7429-90-5]	~0,0026	%		
	Arsenic (As)	[7440-38-2]	0,045	%		
	Bismuth (Bi)	[7440-69-9]	0,038	%		
	Cadmium (Cd)	[7440-43-9]	0,0103	%		
	Copper (Cu)	[7440-50-8]	4,49	%		
	Iron (Fe)	[7439-89-6]	0,0138	%		
	Indium (In)	[7440-74-6]	0,0179	%		
	Nickel (Ni)	[7440-02-0]	1,22	%		
	Phosphorus (P)	[7723-14-0]	~0,002	%		
	Lead (Pb)	[7439-92-1]	0,056	%		
	Antimony (Sb)	[7440-36-0]	4,81	%		
	Selenium (Se)	[7782-49-2]	0,0038	%		
	Zinc (Zn)	[7440-66-6]	~0,02	%		