

Tin-base Lead-free Solder, 50 g Chips

Art. ID MBH-74XHNF-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Copper (Cu)	[7440-50-8]	4,12	%		
	Silver (Ag)	[7440-22-4]	0,16	%		
	Lead (Pb)	[7439-92-1]	0,005	%		
	Antimony (Sb)	[7440-36-0]	0,038	%		
	Bismuth (Bi)	[7440-69-9]	0,042	%		
	Zinc (Zn)	[7440-66-6]	0,0068	%		
	Nickel (Ni)	[7440-02-0]	0,195	%		
	Iron (Fe)	[7439-89-6]	0,02	%		
	Arsenic (As)	[7440-38-2]	0,016	%		
	Cadmium (Cd)	[7440-43-9]	0,0046	%		
	Selenium (Se)	[7782-49-2]	0,0024	%		
	Aluminium (Al)	[7429-90-5]	0,0011	%		
	Indium (In)	[7440-74-6]	0,0052	%		
	Phosphorus (P)	[7723-14-0]	~0,001	%		