

Tin-base (Lead-Free) Solder, 50 g Chips

Art. ID MBH-74XTCF-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,039	%		
	Arsenic (As)	[7440-38-2]	0,024	%		
	Bismuth (Bi)	[7440-69-9]	0,106	%		
	Cadmium (Cd)	[7440-43-9]	0,015	%		
	Copper (Cu)	[7440-50-8]	4,99	%		
	Iron (Fe)	[7439-89-6]	0,0031	%		
	Indium (In)	[7440-74-6]	0,0215	%		
	Nickel (Ni)	[7440-02-0]	0,0167	%		
	Lead (Pb)	[7439-92-1]	0,183	%		
	Antimony (Sb)	[7440-36-0]	0,124	%		
	Selenium (Se)	[7782-49-2]	0,0473	%		
	Zinc (Zn)	[7440-66-6]	0,004	%		