

## **Sn/Pb Solder, 50 g Chips**

Art. ID MBH-91XS40PR2D-C  
Unit 50 g (chips)  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,086	%		
	Arsenic (As)	[7440-38-2]	0,01	%		
	Bismuth (Bi)	[7440-69-9]	0,154	%		
	Cadmium (Cd)	[7440-43-9]	0,0046	%		
	Copper (Cu)	[7440-50-8]	0,085	%		
	Iron (Fe)	[7439-89-6]	0,0096	%		
	Nickel (Ni)	[7440-02-0]	0,005	%		
	Antimony (Sb)	[7440-36-0]	0,596	%		
	Tin (Sn)	[7440-31-5]	40,68	%		
	Zinc (Zn)	[7440-66-6]	0,0275	%		