

Sn/Pb Solder, 50 g Chips

Art. ID MBH-91XS63PR2K-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0528	%		
	Aluminium (Al)	[7429-90-5]	0,0005	%		
	Arsenic (As)	[7440-38-2]	0,0094	%		
	Gold (Au)	[7440-57-5]	0,082	%		
	Bismuth (Bi)	[7440-69-9]	0,175	%		
	Cadmium (Cd)	[7440-43-9]	0,0089	%		
	Copper (Cu)	[7440-50-8]	0,0841	%		
	Iron (Fe)	[7439-89-6]	0,0053	%		
	Indium (In)	[7440-74-6]	0,0157	%		
	Nickel (Ni)	[7440-02-0]	0,0027	%		
	Antimony (Sb)	[7440-36-0]	0,531	%		
	Tin (Sn)	[7440-31-5]	61,84	%		
	Tellurium (Te)	[13494-80-9]	0,0036	%		
	Zinc (Zn)	[7440-66-6]	0,0032	%		