

Sn/Pb/Sb Solder, 50 g Chips

Art. ID MBH-93XS30APR1C-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,0144 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,01 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,059 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,0014 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,192 | % | | |
| | Indium (In) | [7440-74-6] | 0,0094 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,001 | % | | |
| | Antimony (Sb) | [7440-36-0] | 2,54 | % | | |
| | Tin (Sn) | [7440-31-5] | 28,58 | % | | |