

Sn/Pb/Sb Solder, 50 g Chips

Art. ID MBH-93XS30APR2C-C
Unit 50 g
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,049	%		
	Arsenic (As)	[7440-38-2]	0,0178	%		
	Bismuth (Bi)	[7440-69-9]	0,168	%		
	Cadmium (Cd)	[7440-43-9]	0,0061	%		
	Copper (Cu)	[7440-50-8]	0,062	%		
	Iron (Fe)	[7439-89-6]	0,0026	%		
	Indium (In)	[7440-74-6]	0,0199	%		
	Nickel (Ni)	[7440-02-0]	0,042	%		
	Antimony (Sb)	[7440-36-0]	1,8	%		
	Tin (Sn)	[7440-31-5]	30,68	%		
	Tellurium (Te)	[13494-80-9]	0,0102	%		
	Zinc (Zn)	[7440-66-6]	0,028	%		