

Mg/Mn Alloy, chippings 50 g

Art. ID MBH-63XMGE3C-C
Unit 50 g (chips)
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	~0,056	%		
	Zinc (Zn)	[7440-66-6]	0,0101	%		
	Manganese (Mn)	[7439-96-5]	1,62	%		
	Copper (Cu)	[7440-50-8]	0,0072	%		
	Silicon (Si)	[7440-21-3]	0,014	%		
	Iron (Fe)	[7439-89-6]	0,0005	%		
	Nickel (Ni)	[7440-02-0]	0,0028	%		
	Tin (Sn)	[7440-31-5]	0,006	%		
	Lead (Pb)	[7439-92-1]	0,0046	%		
	Silver (Ag)	[7440-22-4]	0,0041	%		
	Beryllium (Be)	[7440-41-7]	<0,0001	%		
	Cadmium (Cd)	[7440-43-9]	0,0046	%		