

Medium Test Dust (MTD) in Hydraulic Fluid (NIST-2806d should only be used in conjunction with ISO method ISO 11171:2020 or a later version) - ONLY LIMITED AVAILABILITY

Art. ID	NIST-2806d
Unit	3 x 400 mL
Deliverydetails	No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM) is intended for use in the calibration of instrumental response to medium test dust suspended in hydraulic fluid. The reference material should only be used with the document ISO 11171:2020 or a later version. A unit of NIST-2806d consists of three bottles containing polydisperse, irregularly shaped mineral dust suspended in approximately 400 mL of hydraulic fluid each.