

**Out of Stock - Product is not available at this time - Copper Freezing-Point Reference (1084.6 Degrees Celsius)**

Art. ID                      NIST-8172  
Unit                         450 g  
Deliverydetails            No Dangerous Good /not restricted

Description

Reference material (RM) 8172 is intended primarily for use as a reference point for allowing laboratories to compare their present realization of a temperature scale at the copper freezing-point with NIST's historical realizations of the International Practical Temperature Scale of 1968 (IPTS-68) circa December 6, 1971 and of the International Temperature Scale of 1990 (ITS-90) circa April 12, 1990, where the dates correspond to the dates on the certificates initially provided with SRM 45d and when updated for the transition from IPTS-68 to ITS-90.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Freezing point		1084,6	Grad Celsi us		