

**Out of Stock - Product is not available at this time - Human Liver for Proteomics**

Art. ID	NIST-8461
Unit	0.5 g
Deliverydetails	No Dangerous Good /not restricted

**Description**

This Reference Material (RM) is intended to support investigations of peptide and protein identification in complex samples such as those in mass spectrometry-based proteomics. RM 8461 can be used to help assess the confidence of peptide and protein identification within a laboratory or comparability between laboratories or among different measurement approaches. RM 8461 can also be used in the development and validation of new investigative approaches for identifying peptides and proteins. A unit of RM 8461 consists of a single vial containing approximately 0.5 g of freeze-dried human liver material.