

Iron Base, Steel - Resulfurized, CRM, Disc 32-35 mm x 19 mm (approx.) - low stock - please ask for availability

Art. ID VASKUT-B12-20MM
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|---------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,007 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,011 | % | | |
| | Boron (B) | [7440-42-8] | 0,0035 | % | | |
| | Carbon (C) | [7440-44-0] | 0,43 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,011 | % | | |
| | Chromium (Cr) | [7440-47-3] | 1,32 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,41 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,76 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,21 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,022 | % | | |
| | Nickel (Ni) | [7440-02-0] | 1,62 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,028 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,048 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,34 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,032 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,011 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,026 | % | | |
| informational value | Zirconium (Zr) | [7440-67-7] | 0,002 | % | | |