

Aluminium base, Cu/Fe/Mg/Si, SUS, Disc 45 mm x 35 mm (approx.)

Art. ID VASKUT-ASC-1-35MM
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
approximate analysis	Beryllium (Be)	[7440-41-7]	0,003	%		
approximate analysis	Bismuth (Bi)	[7440-69-9]	0,1	%		
approximate analysis	Calcium (Ca)	[7440-70-2]	0,02	%		
approximate analysis	Cadmium (Cd)	[7440-43-9]	0,05	%		
approximate analysis	Chromium (Cr)	[7440-47-3]	0,2	%		
approximate analysis	Copper (Cu)	[7440-50-8]	6	%		
approximate analysis	Iron (Fe)	[7439-89-6]	1,6	%		
approximate analysis	Gallium (Ga)	[7440-55-3]	0,04	%		
approximate analysis	Magnesium (Mg)	[7439-95-4]	1,2	%		
approximate analysis	Manganese (Mn)	[7439-96-5]	0,4	%		
approximate analysis	Nickel (Ni)	[7440-02-0]	0,6	%		
approximate analysis	Lead (Pb)	[7439-92-1]	0,1	%		
approximate analysis	Antimony (Sb)	[7440-36-0]	0,02	%		
approximate analysis	Silicon (Si)	[7440-21-3]	14	%		
approximate analysis	Tin (Sn)	[7440-31-5]	0,1	%		
approximate analysis	Strontium (Sr)	[7440-24-6]	0,03	%		
approximate analysis	Titanium (Ti)	[7440-32-6]	0,5	%		
approximate analysis	Vanadium (V)	[7440-62-2]	0,02	%		
approximate analysis	Zinc (Zn)	[7440-66-6]	0,5	%		
approximate analysis	Zirconium (Zr)	[7440-67-7]	0,01	%		