

Custom Solution - 5 components 100 mg/L each of in HNO3 2%, HF tr., C4H6O6 tr.

Art. ID LM24-CUS-71220
Unit 1 x 125 mL bottle
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 III E1

Description

Conservation 8 months (originally closed bottle), (4 months after opening).

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	100	mg/L		
	Scandium (Sc)	[7440-20-2]	100	mg/L		
	Germanium (Ge)	[7440-56-4]	100	mg/L		
	Yttrium (Y)	[7440-65-5]	100	mg/L		
	Bismuth (Bi)	[7440-69-9]	100	mg/L		