

Copper base, Type: Low Alloy, CRM, Disc 40 mm x 30 mm

Art. ID IMN-CM3
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0029	%		
	Arsenic (As)	[7440-38-2]	0,0036	%		
	Bismuth (Bi)	[7440-69-9]	0,0033	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,012	%		
	Nickel (Ni)	[7440-02-0]	0,0031	%		
	Phosphorus (P)	[7723-14-0]	0,0041	%		
	Lead (Pb)	[7439-92-1]	0,0038	%		
	Antimony (Sb)	[7440-36-0]	0,004	%		
	Tin (Sn)	[7440-31-5]	1,06	%		
	Zinc (Zn)	[7440-66-6]	0,006	%		