

**Cylindrospermopsin (CYN), approx. 10 µg/mL in deionized water (ISO 17034)**

Art. ID                      CFG-CRM-03-CYN  
Unit                         0,5 mL  
Deliverydetails            Dry ice shipment /not restricted

Description

EXAMPLE CONCENTRATION: [CYN] = (10.8 ± 0.5) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cylindrospermopsin	[143545-90-8]	approx. 10	µg/mL		