

Gonyautoxins 2 & 3 (GTX2 and GTX3), approx. 30 µg/mL Aqueous [AcOH] = 25·10⁻⁶ M, pH = 4.92 (ISO 17034)

Art. ID CFG-CRM-00-GTX2and3
Unit 0,5 mL
Deliverydetails De Minimis

Description

EXAMPLE CONCENTRATION: [GTX2] = (22.2 ± 1.5) µg/g [GTX3] = (8.2 ± 0.6) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gonyautoxin 2	[60508-89-6]	approx. 30	µg/mL		