

**Okadaic acid (OA), approx. 6 µg/mL in Methanol (ISO 17034)**

Art. ID                      CFG-CRM-00-OA  
Unit                         0,5 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Description

EXAMPLE CONCENTRATION: (20.2 ± 1.0) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Okadaic acid sodium sal t	[209266-80-8]	approx. 16	µg/mL		