

Azaspiracid-3 (AZA3), approx. 1 µg/mL in Methanol (ISO 17034)

Art. ID CFG-CRM-02-AZA3
Unit 0,5 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Description

EXAMPLE CONCENTRATION: AZA3+37-epi-AZA3] = (1.32± 0.09) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Azaspiracid-3	[265996-93-8]	approx. 1	µg/mL		