

**Gonyautoxin 5 (GTX5), approx. 20 µg/mL in aqueous [AcOH] = 25 µM, pH = 5.00 (ISO 17034)**

Art. ID                      CFG-CRM-00-GTX5  
Unit                         0,5 mL  
Deliverydetails            De Minimis

Description

EXAMPLE CONCENTRATION: [GTX5] = (18.1 ± 1.2) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gonyautoxin-5	[64296-25-9]	approx. 20	µg/mL		