

**Gonyautoxin 6 (GTX6), approx. 10 µg/mL in aqueous [AcOH] = 25 µM, pH = 5.00 (ISO 17034)**

Art. ID                      CFG-CRM-00-GTX6  
Unit                         0,5 mL  
Deliverydetails            De Minimis

Description

EXAMPLE CONCENTRATION: [GTX6] = (10.0 ± 0.5) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Gonyautoxin-6	[82810-44-4] (free base)	approx. 10	µg/mL		