

**7-O-Palmitoyl Okadaic Acid (DTX3a), approx. 6 µg/mL in Methanol (ISO 17034)**

Art. ID                      CFG-CRM-02-DTX3a  
Unit                         0,5 mL  
Deliverydetails            De Minimis UN1230 (6.1) II E2

Description

EXAMPLE CONCENTRATION: (7.97 ± 0.48) µg/g

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	7-O-palmitoyl okadaic acid sodium salt	[209266-80-8][57-10-3]	approx. 6	µg/mL		