

**Zearalenone (ZON), 50 µg/mL in Acetonitrile (ISO 17034)**

Art. ID                      CFG-CRM-00-ZON  
Unit                         1 mL  
Deliverydetails            De Minimis UN1648 II E2

Description

EXAMPLE CONCENTRATION: [ZON] = (65 ± 5) µg/g

| Text/Information | Analyte/Parameter | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|--------------|---------------------|-------|--------|--------|
|                  | Zearalenone       | [17924-92-4] | approx. 50          | µg/mL |        |        |