

**Copper base, Type: Phosphorus Alloy, Disc 40 mm x 25 mm**

Art. ID IMN-CO3  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,011	%		
	Bismuth (Bi)	[7440-69-9]	0,015	%		
	Copper (Cu)	[7440-50-8]	Rem	%		
	Iron (Fe)	[7439-89-6]	0,11	%		
	Nickel (Ni)	[7440-02-0]	0,1	%		
	Phosphorus (P)	[7723-14-0]	8,56	%		
	Lead (Pb)	[7439-92-1]	0,1	%		
	Antimony (Sb)	[7440-36-0]	0,14	%		
	Selenium (Se)	[7782-49-2]	~0,0073	%		
	Tin (Sn)	[7440-31-5]	0,037	%		
	Tellurium (Te)	[13494-80-9]	~0,008	%		
	Zinc (Zn)	[7440-66-6]	0,24	%		