

## **Zinc-copper disc, 38x32mm**

Art. ID DHJ-58A-CS04299  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Silver (Ag)       | [7440-22-4] | 0,00096             | %    |        |        |
|                  | Aluminium (Al)    | [7429-90-5] | 0,084               | %    |        |        |
|                  | Arsenic (As)      | [7440-38-2] | 0,044               | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,0044              | %    |        |        |
|                  | Cadmium (Cd)      | [7440-43-9] | 0,00017             | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 62,24               | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,0076              | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,097               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,11                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0077              | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,012               | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,0046              | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 0,03                | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,12                | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | rem                 | %    |        |        |