

Pure copper disc, 40x30mm

Art. ID DHJ-58A-YA04008
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,003 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,00057 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,00014 | % | | |
| | Cadmium (Cd) | [7440-43-9] | 0,00046 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,00056 | % | | |
| | Chromium (Cr) | [7440-47-3] | 0,00046 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,00066 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,00042 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,00061 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,00085 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,00055 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0019 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,00028 | % | | |
| | Selenium (Se) | [7782-49-2] | 0,00019 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,00054 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,00055 | % | | |
| | Tellurium (Te) | [13494-80-9] | 0,00015 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,00093 | % | | |