

Pure copper disc, 40x30mm

Art. ID DHJ-58A-CP04215
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Only available as set. Please see DHJ-58A-CP04213,4,5,6,7

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,042 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,0032 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,0031 | % | | |
| | Cu+Ag | | 99,64 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,046 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,194 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0274 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,0022 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0013 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,0086 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,0318 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,0065 | % | | |