

Pure copper disc, 40x30mm

Art. ID DHJ-58A-CP04217
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Only available as set. Please see DHJ-58A-CP04213,4,5,6,7

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|--------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 0,011 | % | | |
| | Arsenic (As) | [7440-38-2] | 0,0158 | % | | |
| | Bismuth (Bi) | [7440-69-9] | 0,0069 | % | | |
| | Cadmium (Cd) | [7440-43-9] | ~0,0075 | % | | |
| | Cu+Ag | | 99,82 | % | | |
| | Iron (Fe) | [7439-89-6] | 0,0109 | % | | |
| | Manganese (Mn) | [7439-96-5] | ~0,0007 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,0158 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0014 | % | | |
| | Lead (Pb) | [7439-92-1] | 0,0143 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0151 | % | | |
| | Antimony (Sb) | [7440-36-0] | 0,016 | % | | |
| | Selenium (Se) | [7782-49-2] | ~0,0074 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,015 | % | | |
| | Tellurium (Te) | [13494-80-9] | ~0,0079 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,002 | % | | |