

Copper base, Type: Brass - Manganese, CRM, Set IMN-MA1 - IMN-MA6 only, see IMN-MA1-MA6, Rod 10 mm x 100 mm

Art. ID IMN-MA3
Unit rod
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	0,029	%		
	Bismuth (Bi)	[7440-69-9]	0,028	%		
	Copper (Cu)	[7440-50-8]	57,04	%		
	Iron (Fe)	[7439-89-6]	0,55	%		
	Manganese (Mn)	[7439-96-5]	0,78	%		
	Nickel (Ni)	[7440-02-0]	0,13	%		
	Phosphorus (P)	[7723-14-0]	0,04	%		
	Lead (Pb)	[7439-92-1]	0,049	%		
	Antimony (Sb)	[7440-36-0]	0,14	%		
	Silicon (Si)	[7440-21-3]	0,5	%		
	Tin (Sn)	[7440-31-5]	0,74	%		
	Zinc (Zn)	[7440-66-6]	rem	%		