

**Copper base, Type: Brass, CRM, wrought, Disc 40 mm x 30 mm**

Art. ID IMN-MG6  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Aluminium (Al)    | [7429-90-5] | 0,0067              | %    |        |        |
|                  | Bismuth (Bi)      | [7440-69-9] | 0,00088             | %    |        |        |
|                  | Copper (Cu)       | [7440-50-8] | 92,27               | %    |        |        |
|                  | Iron (Fe)         | [7439-89-6] | 0,028               | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,045               | %    |        |        |
|                  | Nickel (Ni)       | [7440-02-0] | 0,03                | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,0026              | %    |        |        |
|                  | Lead (Pb)         | [7439-92-1] | 0,031               | %    |        |        |
|                  | Antimony (Sb)     | [7440-36-0] | 0,0015              | %    |        |        |
|                  | Tin (Sn)          | [7440-31-5] | 0,053               | %    |        |        |
|                  | Zinc (Zn)         | [7440-66-6] | rem                 | %    |        |        |