

**Copper base, Type: Brass - Cartridge, CRM, wrought, Disc 40 mm x 25 mm**

Art. ID IMN-MH1  
Unit disc  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,0029	%		
	Aluminium (Al)	[7429-90-5]	0,001	%		
	Arsenic (As)	[7440-38-2]	0,067	%		
	Beryllium (Be)	[7440-41-7]	0,0088	%		
	Bismuth (Bi)	[7440-69-9]	0,0037	%		
	Cadmium (Cd)	[7440-43-9]	0,026	%		
	Copper (Cu)	[7440-50-8]	65,93	%		
	Iron (Fe)	[7439-89-6]	0,017	%		
	Manganese (Mn)	[7439-96-5]	0,035	%		
	Nickel (Ni)	[7440-02-0]	0,26	%		
	Phosphorus (P)	[7723-14-0]	0,016	%		
	Lead (Pb)	[7439-92-1]	0,0065	%		
	Sulfur (S)	[7704-34-9]	0,0034	%		
	Antimony (Sb)	[7440-36-0]	0,0004	%		
	Silicon (Si)	[7440-21-3]	0,074	%		
	Tin (Sn)	[7440-31-5]	0,14	%		
	Tellurium (Te)	[13494-80-9]	0,0004	%		
	Zinc (Zn)	[7440-66-6]	rem	%		