

Copper base, Type: Gilding Metal, CRM, wrought, Disc 40 mm x 25 mm

Art. ID IMN-MI5
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	0,033	%		
	Aluminium (Al)	[7429-90-5]	0,0021	%		
	Arsenic (As)	[7440-38-2]	0,015	%		
	Beryllium (Be)	[7440-41-7]	0,0072	%		
	Bismuth (Bi)	[7440-69-9]	0,0043	%		
	Cadmium (Cd)	[7440-43-9]	0,0012	%		
	Copper (Cu)	[7440-50-8]	94,71	%		
	Iron (Fe)	[7439-89-6]	0,015	%		
	Manganese (Mn)	[7439-96-5]	0,069	%		
	Nickel (Ni)	[7440-02-0]	0,25	%		
	Phosphorus (P)	[7723-14-0]	0,0026	%		
	Lead (Pb)	[7439-92-1]	0,096	%		
	Sulfur (S)	[7704-34-9]	0,0019	%		
	Antimony (Sb)	[7440-36-0]	0,0096	%		
	Silicon (Si)	[7440-21-3]	0,082	%		
	Tin (Sn)	[7440-31-5]	0,004	%		
	Zinc (Zn)	[7440-66-6]	4,44	%		